



The biotIP® Confinement Mask, SEP4 for use with multiple sensing electrode platform, SEP1.

The biotIP® confinement maskSEP4, is used in conjunction with the multichannel sensing test platform SEP1.

Made from high quality hydrophobic substrate, the self-adhesive mask is designed specifically to isolate any individual channel for more accurate and easier drop casting.

The self-adhesive masks are laser cut and precision engineered to fit directly over each test channel for parallel drop casting as demonstrated below.




Figure 1: biotIP® SEP4, confinement mask.



Effortless and accurate alignment is accomplished by the use of laser cut nodes which match on both the test platform as well as the mask for consistent confinement and coverage.




Figure 2: biotIP® SEP4, confinement mask fitting over the SEP1 electrode set.